

Introducing Surestep Go SMO

Surestep Go SMO is a prefabricated SMO designed to manage mild to moderate pronation associated with hypotonia. This SMO uses circumferential compression to stabilize the foot and ankle. The open heel allows more of the foot to "feel the floor", which aids in sensory input and posture.

For maximum function, this SMO should always be worn with footwear



Support for Better Life!

PRODUCT FEATURES

Total contact/circumferential compression

Surestep Go SMO provides circumferential compression to manage excessive pronation and provide dynamic stability.

Surestep Trimlines

The trimlines on the footplate are posterior to the forefoot (metatarsal heads) on the medial side. This frees the toes to push off and allows them to feel the floor.

Open Heel

The open heel improves proprioception by allowing the heel to feel the floor and encourages posterior weight shifts.

AFO & SMO Combination

Surestep Go SMO can be combined with KiddieGAIT®, KiddieROCKER® and KiddieFLOW™ to support gait deviations and manage Hypotonia with mild to moderate foot pronation.



INDICATIONS

Mild to moderate pronation, hypotonia, triplanar instability in weight bearing, inability to stand independently, mild toe walking

CONTRAINDICATIONS

Hypertonia, Severe Spasticity, Severe Toe Walking

Item No. Pair	Item No. Left	Item No. Right	Size	UK Shoe Size	EU Shoe Size
28740 0003	28742 1003	28742 2003	3T	2	19
28740 0004	28742 1004	28742 2004	4T	4	20
28740 0005	28742 1005	28742 2005	5T	5	21
28740 0006	28742 1006	28742 2006	6T	5.5	22
28740 0007	28742 1007	28742 2007	7T	6	23
28740 0008	28742 1008	28742 2008	8T	7	24
28740 0009	28742 1009	28742 2009	9T	8.5	26
28741 0010	28743 1010	28743 2010	10Y	9	27
28741 0011	28743 1011	28743 2011	11Y	10	28
28741 0012	28743 1012	28743 2012	12Y	11	29

Support for Better Life!