

Elements Body Sock

Item no: 39325

A custom made dynamic textile Ankle Foot Orthosis. The Elements Body Sock is a closed toe, below knee design, with an anterior zip fastening. This garment is designed to assist movement of the foot and ankle complex, whilst providing stabilisation, and sensory and proprioceptive feedback. This is achieved by adding strategic reinforcements that provide correction and directional pull.

Recommended Clinical Application:

Elements Body Sock may be prescribed for the following conditions where there is passive range of movement in the ankle and foot and some activity against gravity:

- Cerebral Palsy
- Acquired Brain Injury
- Acquired Spinal Cord Injury
- Cerebellar Ataxia
- Athetosis
- Spina Bifida
- Stroke
- Multiple Sclerosis
- Hyper-mobility

The Elements Body Sock may also be used in the management of joint instability and/or recovering peripheral nerve injury.

Contraindications:

- Lower Limb Ulcers
- Vascular Disorders

