

INFORMATION ABOUT SOME MEASUREMENTS

If the neck line is being made deeper or wider please take note of the following maximum extremities:

For children, measurement ND cannot exceed 3cm. For adults, ND cannot exceed 6cm. For children, measurement NS cannot exceed 1cm. For adults, NS cannot exceed 3cm. For children, measurement BD cannot exceed 3cm. For adults, BD cannot exceed 4cm.

It is now possible to include a Velcro tab to attach a Vest to Trousers

When ordering, the measurement UE needs to be extended at least 3.5cm below ASIS on the Vest for Children. For Adults, the UE measurement should be extended at least 4.5cm below ASIS on the Vest.