

allard_{UK}

ToeOFF[®] SHORT

BLUEROCKER[®] SHORT

New Model!



NEW features

UPDATED STRAP SOLUTION
& PRE-APPLIED MICROFIX!

ToeOFF[®] SHORT BLUEROCKER[®] SHORT

SHORT models are specifically developed to accommodate individuals with a shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area. Tibia padding should always be worn for patient comfort.

Both SHORT models have a heel height of 7 mm in all sizes.

Intended Use

Foot drop in combination with weak plantar flexors, and in some cases limited knee control.

For ToeOFF[®] SHORT: Conditions such as Stroke, Multiple Sclerosis, Post Polio Syndrome, Muscular Dystrophy, and other neurological diagnosis.

For BlueROCKER[®] SHORT: This may also be used for bilateral problems and partial foot amputations, most proximal level is Chopart.

Contraindications

Should not be used when patients present with foot and/or leg ulcers, moderate to severe edema, moderate to severe foot deformities, severe proximal deficits or severe spasticity.

More about fitting, indications and product selection can be found in the Allard AFO Professional Instructions.

Item No. ToeOFF [®]	Item No. BlueROCKER [®]	Size	Side	Total height	Footplate length
28990 1011	28995 1011	Small	Left	30,5 cm	21,5 cm
28990 1012	28995 1012	Medium	Left	32,5 cm	23,2 cm
28990 1013	28995 1013	Large	Left	34 cm	25 cm
28990 2011	28995 2011	Small	Right	30,5 cm	21,5 cm
28990 2012	28995 2012	Medium	Right	32,5 cm	23,2 cm
28990 2013	28995 2013	Large	Right	34 cm	25 cm
28391 0011	SoftKIT [™] SHORT, Small				
28391 0012	SoftKIT [™] SHORT, Medium				
28391 0013	SoftKIT [™] SHORT, Large				
28768 0011	SHORT Strap set, Small				
28768 0012	SHORT Strap set, Medium				
28768 0013	SHORT Strap set, Large				

