

SELECTION[®]

CHILDREN'S OPEN WRIST



Wrist Orthosis for Children

Intended to help wrist and adducted thumb to a neutral and more functional mode.

- Thumb in a neutral position
- Support for the wrist
- Soft and comfortable material
- Extra stabilizing strap around the wrist.

SELECTION[®]

CHILDREN'S OPEN WRIST



The two-piece thumb strap keeps the thumb in a neutral and functional position.



Integrated stabilizing strap around the wrist.



The extra material supports the thumb's muscles and the extra strap helps to move the adducted thumb out from the palm.

Item No. left	Item No. right	Size	Wrist circumference	Palm length
35306 1011	35306 2011	Small	11,5 - 14 cm	7 - 8 cm
35306 1012	35306 2012	Medium	12 - 15,5 cm	8 - 9 cm
35306 1013	35306 2013	Large	14 - 16,5 cm	8 - 9 cm

