

Stock Size Chart

Boston Scoliosis, Soft Brace, Boston LITE, Flexaform, B.O.B, Soft B.O.B.

Size chart B.O.B

Size	Hip cm	Waist cm	Chest cm
1	71	51	66
1D	71	45	58,5
2	76	45	58,5
2D	76	51	63,5
3	81	71	76
3D	81	66	74
4	86,5	76	84
4D	86,5	58,5	69
5	91,5	81	89
5D	91,5	66	76
6	96,5	86,5	94
6D	96,5	74	84
7	101,5	91,5	99
7D	101,5	81	91,5
8	109	96,5	109
8D	109	84	91,5
9	114	100	114
9D	114	87	94
10	117	104	117
10D	117	96	101

Size chart Scoliosis and Soft Body Jacket

Size	Hip cm	Waist cm	Chest cm
J1	54	40	49
J2	58	42	51
J3	62	45	54
J4	64	47	56
Y1A	66	51	55
Y1	67	46	56
Y2	71	47	60
Y2A	72	55	62
Y3	75	48	63
Y3A	78	57	69
Y4	79	51	64
Y4F	80	57	66
Y5	82	52	68
Y5F	82	58	68
Y6	84	56	73
Y6A	84	63	78
Y7	86	54	69
Y7F	86	60	69
Y7A	86	69	73
Y8	88	58	75
Y8F	88	64	75
Y8A	90	72	83
Y9	92	60	72
Y9A	92	63	69
Y9F	92	68	78
A1	94	62	76
A2	94	70	84
A2A	96	66	75
A3	97	68	80
A4	97	76	88

Size chart Flexaform

Size	Hip cm	Waist cm	Chest cm
LS1	79	60	71
LS2	82	65	76
LS3	85	63	75
LS4	87	69	82
LS5	90	68	80
LS6	92	75	88
LS7	95	74	85
LS8	97	83	94
LS9	100	80	90
LS10	102	87	99
LS11	105	87	95
LS12	107	93	104

Size Chart Boston LITE

Item No.	Size	Hip cm	Waist cm	Chest cm
SBBL-15J1	15J1	54	40	49
SBBL-15J2	15J2	58	42	51
SBBL-15J3	15J3	62	45	54
SBBL-15J4	15J4	64	47	56
SBBL-15Y1	15Y1	67	46	56
SBBL-15Y1A	15Y1A	66	51	55
SBBL-15Y2	15Y2	71	47	60
SBBL-15Y2A	15Y2A	72	55	62
SBBL-15Y3	15Y3	75	48	63
SBBL-15Y3A	15Y3A	78	57	69
SBBL-15Y4	15Y4	79	51	64
SBBL-15Y4F	15Y4F	80	57	66

Please quote the size and lordosis (0° or 15°) when ordering.